

## PTA BALLOONS

| Company Name                   | Product Name                                      | Nominal (atm)            | Rated Burst (atm)                | Endhole (inch)                   | Shaft Length (cm)            | Sheath Sizing Size (F)                              | Balloon Diameters (mm)           | Balloon Lengths (cm)                |
|--------------------------------|---|--------------------------|----------------------------------|----------------------------------|------------------------------|---|----------------------------------|-------------------------------------|
| Abbott Vascular                | Armada 14 PTA Catheter                            | 8                        | 14                               | N/A                              | Catheter length: 90, 150     | Minimum sheath size: 4                              | 1.5, 2, 2.5, 3, 4                | 20, 40, 60, 80, 120, 200 mm         |
|                                | Armada 35 LL PTA Catheter                         | 6, 8                     | 12, 13, 14, 16, 18, 20           |                                  | Catheter length: 80, 135     | Minimum sheath size: 5, 6                           | 4, 5, 6, 7                       | 150, 200, 250 mm                    |
|                                | RX Viatrac 14 Plus Peripheral Dilatation Catheter | 8                        | 14                               |                                  | Catheter length: 80, 135     | Maximum guidewire: 4, 5                             | 4, 4.5, 5, 5.5, 6, 6.5, 7        | 15, 20, 30, 40 mm                   |
|                                | Agiltrac 0.018 Peripheral Dilatation Catheter     |                          |                                  |                                  |                              | N/A   | 4, 5, 6, 7, 8, 9, 10, 12, 14     | 20, 40, 60, 80, 100 mm              |
|                                | Fox Plus  | 5, 7, 8, 9, 10           | 8, 9, 12, 13, 14, 15, 16, 17, 18 |                                  | Catheter length: 50, 80, 135 | Maximum guidewire: 0.035-inch                       | 3, 4, 5, 6, 7, 8, 9, 10, 12      | 20, 30, 40, 60, 80, 100, 120 mm     |
|                                | Fox sv  | 6, 8, 12                 | 15, 16, 17, 18, 20, 22           |                                  | Catheter length: 90, 150     | Maximum guidewire: 0.018-inch                       | 2, 2.5, 3, 4, 5, 6               | 15, 20, 30, 40, 60, 80, 100, 120 mm |
|                                | FoxCross  | N/A                      | Up to 18                         |                                  | Catheter length: 50, 80, 135 | Maximum guidewire: 0.035-inch                       | 3, 4, 5, 6, 7, 8, 9, 10, 12, 14  | 20, 30, 40, 60, 80, 100, 120 mm     |
| AngioDynamics                  | WorkHorse II                                      | 10                       | 14, 16, 18                       | 0.035                            | 50, 75, 100                  | 5–8 mm through a 6 F, 9–10 mm through a 7 F         | 5, 6, 7, 8, 9, 10                | 2, 4, 6                             |
|                                | Profiler  | 7                        | 12, 14, 16                       |                                  | 40, 75, 120                  | 4–8 X 4 mm through a 5 F, 8 X 6–10 mm through a 6 F | 4, 5, 6, 7, 8, 9, 10             | 2, 4, 6, 8                          |
| AngioScore Inc.                | AngioSculpt                                       | 8                        | Up to 20                         | 0.014                            | 137                          | 5   | 2, 2.5, 3, 3.5                   | 1, (2-mm only), 2, 4                |
|                                |   |                          | Up to 16                         | 0.018                            | 50 (6 mm only), 90, 137      | 6   | 4, 5, 6                          | 2, 4                                |
|                                | AngioSculpt XL                                    | 14                       | 0.014                            | 90, 137                          |                              |   | 10                               |                                     |
| Angioslide, Inc.               | Proteus PTA Catheter With Embolic Capture Feature | 8                        | 12                               | 0.035                            | 135                          | 4–5-mm through a 6 F, 6-mm through a 7 F            | 4, 5, 6                          | 6, 8, 10                            |
| Bard Peripheral Vascular, Inc. | Conquest  | 8                        | Up to 30                         | 0.035                            | 50, 75, 120                  | 6–8   | 5, 6, 7, 8, 9, 10, 12            | 2, 3, 4, 6, 8                       |
|                                | Atlas   | 4, 6                     | Up to 18                         |                                  | 75, 120                      | 7–12  | 12, 14, 16, 18, 20, 22, 24, 26   | 2, 4, 6                             |
|                                | Rival   | 6, 8                     | Up to 16                         |                                  | 80, 135                      | 5–7   | 3, 4, 5, 6, 7, 8, 9, 10          | 2, 4, 6, 8, 10, 15                  |
|                                | Dorado  | 8                        | Up to 24                         |                                  | 40, 80, 120, 135             |   | 3, 4, 5, 6, 7, 8, 9, 10          | 2, 4, 6, 8, 10, 12, 15, 17, 20      |
|                                | Ultraverse 014                                    | 6                        | Up to 16                         | 0.014                            | 150                          | 4–5   | 1.5, 2, 2.5, 3, 3.5, 4, 5        | 2, 4, 8, 12, 15, 22                 |
|                                | Ultraverse 018                                    |                          |                                  | 0.018                            | 75, 130                      | 4–6   | 2, 2.5, 3, 3.5, 4, 5, 6, 7, 8, 9 | 2, 4, 6, 8, 10, 12, 15, 22          |
|                                | VascuTrak   | 8 (2–3.5 mm), 6 (4–7 mm) | 12                               | 0.014 (2–3.5 mm), 0.018 (4–7 mm) | 80, 140                      | 5 (2–4 mm), 6 (5 mm), 7 (6–7 mm)                    | 2, 2.5, 3, 3.5, 4, 5, 6, 7       | 2, 4, 6, 8, 10, 12, 15, 20, 25, 30  |

## PTA BALLOONS (CONTINUED)

| Company Name                    | Product Name                                       | Nominal (atm)                   | Rated Burst (atm)                          | Endhole (inch) | Shaft Length (cm)     | Sheath Sizing Size (F)                              | Balloon Diameters (mm)           | Balloon Lengths (cm)                 |
|---------------------------------|--|---------------------------------|--|----------------|-----------------------|---|----------------------------------|--------------------------------------|
| B. Braun Interventional Systems | Mini Ghost   | Up to 10                        | Up to 20                                   | 0.018          | 40, 80, 120, 150      | 4   | 2–6                              | 1, 2, 3, 4, 6, 8, 10                 |
|                                 | Z-Med  | Up to 6                         | Up to 15                                   | 0.035          | 80, 100, 120          | 5–16  | 2–40                             | 2, 3, 4, 5, 6                        |
|                                 | Impact   |                                 | Up to 10                                   |                | 75, 120               | 8–14  | 12–25                            | 2, 4                                 |
|                                 | CoEfficient  |                                 | Up to 13                                   | 0.018          | 75, 135               | 4–7   | 4–12                             | 2, 3, 4                              |
| Boston Scientific Corporation   | Mustang Balloon Catheter                           | 10 (3–8 mm, 12 mm), 8 (9–10 mm) | Up to 24                                   | 0.035          | 40, 75, 135           | 5 (3–7 mm), 6 (8–10 mm), 7 (12 mm)                  | 3–10, 12                         | 2, 3, 4, 6, 8, 10, 12, 15, 18, 20    |
|                                 | XXL Balloon Catheter                               | N/A                             | 8 (12, 14 mm), 5 (16, 18 mm)               |                | 75, 120               | 7 (12, 14 mm), N/A (16, 18 mm)                      | 12, 14, 16, 18                   | 2, 4, 6                              |
|                                 | Blue Max Balloon Catheter                          | 6                               | 20 (4–8 mm X 4 cm), 17 (8 mm X 8 cm–10 mm) |                | 40, 75, 120           | 6 (4–6 mm), 7 (7–8 mm), 8/9 (9–10 mm)               | 4, 5, 6, 7, 8, 9, 10             | 2, 3, 4, 8, 10                       |
|                                 | Ultra-Thin Diamond Balloon Catheter                |                                 | 15 (3–8 mm), 12 (9–10 mm), 10 (12 mm)      |                | 40, 75, 120, 135, 150 | 5 (3–6 mm), 6 (7–8 mm), 7 (9–10 mm), 8 (12 mm)      | 3, 4, 5, 6, 7, 8, 9, 10, 12      | 1.5, 2, 3, 4, 6, 8, 10               |
|                                 | Ultra-Thin SDS Balloon Catheter                    | 6–8                             | 12   |                | 50, 75, 90, 135       | 5 (up to 7 X 4), 6 (up to 8 X 4), 7 (9–10 mm)       | 4, 5, 6, 7, 8, 9, 10             | 1.5, 2, 3, 4, 6, 8                   |
|                                 | Synergy Balloon Catheter                           | 6 (3–10 mm), 4 (12 mm)          | 12 (3–8 mm), 10 (9–10 mm), 6 (12 mm)       |                | 50, 75, 90, 135, 150  | 5 (up to 6 X 4), 6 (up to 8 X 4), 7 (9–12 mm)       | 3, 4, 5, 6, 7, 8, 9, 10, 12      | 1.5, 2, 3, 4, 6, 8, 10               |
|                                 | Sterling Monorail PTA Balloon Dilatation Catheter  | 6                               | 14 (3–7 mm), 12 (8 mm)                     |                | 0.018                 | 80, 135   | 4 (up to 8 X 4), 5 (up to 8 X 6) | 3, 3.5, 4, 4.5, 5, 5.5, 6, 6.5, 7, 8 |
|                                 | Sterling Over-the-Wire Balloon Dilatation Catheter |                                 | 14 (4–7 mm), 12 (8 mm), 10 (9–10 mm)       | 40, 80, 135    |                       | 4 (up to 8 X 4), 5 (up to 10 X 4), 6 (up to 10 X 8) | 4, 5, 6, 7, 8, 9, 10             | 2, 3, 4, 6, 8, 10                    |

## PTA BALLOONS (CONTINUED)

| Company Name                              | Product Name  | Nominal (atm) | Rated Burst (atm)            | Endhole (inch) | Shaft Length (cm) | Sheath Sizing Size (F)                            | Balloon Diameters (mm)                 | Balloon Lengths (cm)  |
|---|---|---------------|------------------------------|----------------|-------------------|---|--|---|
| Boston Scientific Corporation (Continued) | Sterling SL Monorail PTA Balloon Dilatation Catheter                  | 6             | 14                           | 0.018          | 90, 150           | 4   | 2, 2.5, 3, 3.5, 4                      | 8, 10, 12, 15   |
|   | Sterling SL Over-the-Wire PTA Balloon Dilatation Catheter             |               |                              |                |                   |   |  | 8, 10, 12, 15 (4-mm diameter on 90-cm shaft not available in 8- or 10-cm balloon lengths) |
|   | Symmetry Balloon Catheter With Glidex Hydrophilic Coating             | 6             | 15                           | 0.018          | 90, 135           | 4 (1.5–4 mm X 4 cm), 5 (4 mm X 10 cm–6 mm X 4 cm) | 1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6 | 2, 4, 10  |
|   | Symmetry Stiff Shaft With Glidex Hydrophilic Coating Balloon Catheter |               |                              |                | 90, 135, 150      |   |  |   |
| Coyote Balloon Catheter (OTW)             | Coyote Balloon Catheter (OTW)   | 8             | 14                           | 0.014          | 90, 150           | 4   | 1.5, 2, 2.5, 3, 3.5, 4                 | 4 (on 1.5 mm only); 6, 8, 10, 12, 15, 22  |
|   | Coyote Balloon Catheter (Monorail)                                    |               |                              |                |                   |   |  |   |
|   | Coyote ES Monorail PTA Balloon Dilatation Catheter                    | 6             | 14 (1.5–3 mm), 12 (3.5–4 mm) | 0.014          | 142–146           | 4   | 1.5, 2, 2.5, 3, 3.5, 4                 | 2, 3, 4   |
|   | Coyote ES Over-the-Wire PTA Balloon Dilatation Catheter               |               |                              |                |                   |   |  |   |
| PolarCath Peripheral Dilatation System    | Peripheral Cutting Balloon (OTW)                                      | 6             | 10                           | 0.018          | 50, 90, 135       | 6 (5, 6 mm), 7 (7, 8 mm)                          | 5–8                                    | 2   |
|   |   |               |                              | 0.014          | 140               | 6   | 2, 2.5, 3, 3.5, 4                      | 1.5   |
|   | Peripheral Cutting Balloon (Monorail)                                 | 6             | 12                           | 0.014          | 140               | 6   | 2, 2.5, 3, 3.5, 4                      | 1.5   |
|   |   |               |                              |                |                   |   |  |   |

## PTA BALLOONS (CONTINUED)

| Company Name                         | Product Name                                  | Nominal (atm)                     | Rated Burst (atm)                   | Endhole (inch) | Shaft Length (cm)    | Sheath Sizing Size (F)                                  | Balloon Diameters (mm)            | Balloon Lengths (cm)               |
|--------------------------------------|---|-----------------------------------|-------------------------------------|----------------|----------------------|---|-----------------------------------|------------------------------------|
| Cook Medical                         | Advance 14LP Low Profile PTA Balloon Catheter | 8                                 | 16                                  | 0.014          | 110, 170             | 4   | 2, 2.5, 3, 4                      | 2, 4, 6, 8, 12, 16, 20             |
|                                      | Advance 18LP Low Profile PTA Balloon Catheter |                                   | 14                                  | 0.018          |                      |   |                                   |                                    |
|                                      | Advance 35LP Low Profile PTA Balloon Catheter | 5-10                              | Up to 15                            | 0.035          | 40, 80, 120          | 5 (up to 8 X 8 mm), 6 (up to 10 X 8 mm), 7 (12 x 14 mm) | 3, 4, 5, 6, 7, 8, 9, 10, 12       | 2, 3, 4, 6, 8, 10, 12, 14          |
|                                      | ATB Advance                                   | 5                                 |                                     |                |                      | 15  | 5-7                               | 4, 5, 6, 7, 8, 9, 10, 12, 14       |
| Cordis Corporation                   | Sleek OTW PTA Dilatation Catheter             | 6                                 | 13-16                               | 0.014          | 150                  | 4   | 1.25-5                            | 1.5-22                             |
|                                      | Savvy Long PTA Dilatation Catheter            |                                   | 15, 13, 12                          | 0.018          | 120, 150             | 4-5   | 2, 2.5, 3, 3.5, 4, 5, 6           | 12, 15, 22                         |
|                                      | Sleek RX PTA Dilatation Catheter              |                                   | 15, 16                              | 0.014          | 150, 155             | 4   | 2, 2.5, 3, 3.5, 4                 | 4, 8, 10, 12, 15, 22               |
|                                      | Opta Pro PTA Dilatation Catheter              | 6, 4 (12 mm)                      | 10, 6 (12 mm)                       | 0.035          | 80, 110, 135         | 5-7   | 3-10, 12                          | 1, 1.5, 2-4, 6, 8, 10              |
|                                      | Powerflex P3 PTA Dilatation Catheter          | 10, 6 (12 mm)                     | 15, 8 (12 mm)                       |                | 40, 65, 80, 110, 135 | 5-8   | 4-10, 12                          | 1, 1.5, 2-4, 6, 8, 10              |
|                                      | Powerflex Extreme PTA Dilatation Catheter     | 15, 12 (10 mm)                    | 20, 17 (9-10 mm)                    |                | 40, 80, 120          | 5-9   | 4-10                              | 2, 4, 6                            |
|                                      | Maxi LD PTA Dilatation Catheter               | 4, 3.5 (22, 25 mm)                | 6, 5 (22, 25 mm)                    |                | 80, 110              | 8-11  | Iliac: 14, 15                     | 2, 4, 6, 8                         |
|                                      | Slalom PTA Dilatation Catheter                | 10 (3-6 mm), 8 (6-7 mm), 6 (8 mm) | 14 (3-5 mm), 12 (6-7 mm), 10 (8 mm) | 0.018          | 80, 120, 135         | 4-6   | 3-8                               | 2, 4                               |
|                                      | Savvy PTA Dilatation Catheter                 | 8 (2-4 mm), 6 (4.5-6 mm)          | 10                                  |                | 80, 120, 150         | 4-5   | 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6 | 2-4, 6, 10                         |
| Aviator Plus PTA Dilatation Catheter | 10  | 14                                | 0.014                               | 142            | 4 (4-6 mm), 5 (7 mm) | 4, 4.5, 5, 5.5, 6, 7                                    | 1.5, 2, 3, 4                      |                                    |
| Covidien/ev3                         | EverCross 0.035 PTA Balloon                   | 10, 7                             | 20, 18, 14, 12, 11, 10              | 0.035          | 40, 80, 135          | 5, 6, 7   | 3, 4, 5, 6, 7, 8, 9, 10, 12       | 1.5, 2, 3, 4, 6, 8, 10, 12, 15, 20 |
|                                      | PowerCross 0.018 PTA Balloon                  | 8                                 | 14                                  | 0.018          | 90, 150              | 4, 5, 6   | 2, 2.5, 3, 4, 5, 6                | 2, 4, 6, 8, 10, 12, 15, 20         |
|                                      | NanoCross 0.014 PTA Balloon                   | 10, 8, 7                          |                                     | 0.014          |                      | 4   | 1.5, 2, 2.5, 3, 3.5, 4            | 2, 4, 8, 12, 15, 21                |

## PTA BALLOONS (CONTINUED)

| Company Name                           | Product Name                            | Nominal (atm) | Rated Burst (atm)          | Endhole (inch) | Shaft Length (cm)         | Sheath Sizing Size (F)   | Balloon Diameters (mm)               | Balloon Lengths (cm)                |
|--|---|---------------|----------------------------|----------------|---------------------------|--------------------------|--------------------------------------|-------------------------------------|
| Hotspur Technologies, Inc.             | IQCath Balloon Dilatation Catheter      | 8             | 20, 16                     | 0.035          | 55                        | 6, 7                     | 5, 6, 7, 8                           | 4                                   |
|  | GPSCath Balloon Dilatation Catheter     |               |                            |                |                           |                          |                                      |                                     |
| Idev (manufactured by TriReme Medical) | Gliderflex                              | 8, 9          | 14, 12                     | 0.014, 0.018   | Catheter length: 120, 135 | 5                        | 2, 2.5, 3, 3.5, 4, 5, 6, 7           | 40, 60, 80, 100, 120, 150, 200, 300 |
| Medtronic, Inc.                        | Amphirion Deep OTW PTA Balloon Catheter | 7             | 14                         | 0.014          | 120, 150                  | 4                        | 1.5, 2, 2.5, 3, 3.5, 4               | 2, 4, 8, 12, 15, 21                 |
|  | Admiral Xtreme OTW PTA Balloon Catheter | 8, 6          | Up to 14                   | 0.035          | 80, 130                   | 5-7                      | 3, 4, 5, 6, 7, 8, 9, 10, 12          | 2, 4, 6, 8, 12, 15, 20, 25, 30      |
|  | Pacific Xtreme PTA Balloon Catheter     | 6             | 12, 14                     | 0.018          | 90, 130                   | 4, 5                     | 4, 5, 6, 7                           | 15, 20, 25, 30                      |
|  | Reef HP OTW PTA Balloon Catheter        | 10            | 22, 20, 18                 |                | 80, 120                   | 5-7                      | 4, 5, 6, 7, 8                        | 2, 3, 4, 6, 8                       |
|  | Submarine Plus OTW PTA Balloon Catheter | 7             | Up to 14                   | 0.035          | 90, 130                   | 4-5                      | 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6, 7 | 2, 4, 6, 8, 12                      |
|  | Optimus PTA Balloon Catheter            | 8             | 12 (5-8 mm), 10 (9, 10 mm) |                | 80, 135                   | 5 (5-8 mm), 6 (9, 10 mm) | 5, 6, 7, 8, 9, 10                    | 20, 30, 40, 60                      |